

MASKEARING

USER MANUAL


MaskEaring is a device created to protect ears from skin irritations and injuries caused by the elastics of face masks if worn for long periods of time.

MaskEaring's aim is to reduce, as much as possible, the discomfort and pain caused by the elastic band of the face mask and, because of that, the attaching points for the elastics are moved forward, to the front of the ears, totally avoiding touching the skin.

MaskEaring is really easy to use and to maintain, just follow the simple steps in the following pages.


HOW TO USE

1. Wear MaskEaring behind your ears, as shown in the image



2. Put the face mask's elastic bands behind the hooks. You can either tie a knot to shorten the elastic bands or you can cross them, in order to fit the mask correctly.

In case you wear glasses or hearing aid, place the temples or the dome between the flaps on the back of the structure of the device. They're really flexible and can easily hold them.



3. Now go on with your day without worrying about your face mask!

HOW TO CLEAN

Whenever you want, you can clean your MaskEaring with a cloth and some detergent or sanitizer.

